

Laboratory Exercise

4

Care and Use of the Microscope

Materials Needed

Compound light microscope
Lens paper
Microscope slides
Coverslips
Transparent plastic millimeter ruler
Prepared slide of letter *e*
Slide of three colored threads
Medicine dropper
Dissecting needle (needle probe)
Specimen examples for wet mounts
Methylene blue (dilute) or iodine-potassium-iodide stain

For Demonstrations:

Micrometer scale
Stereomicroscope (dissecting microscope)

The human eye cannot perceive objects less than 0.1 mm in diameter, so a microscope is an essential tool for the study of small structures such as cells. The microscope usually used for this purpose is the *compound light microscope*. It is called compound because it uses two sets of lenses: an eyepiece or ocular lens system and an objective lens system. The eyepiece lens system magnifies, or compounds, the image reaching it after the image is magnified by the objective lens system. Such an instrument can magnify images of small objects up to about one thousand times.

Purpose of the Exercise

To become familiar with the major parts of a compound light microscope and their functions and to make use of the microscope to observe small objects.

LEARNING OUTCOMES

After completing this exercise, you should be able to

1. Locate and identify the major parts of a compound light microscope and differentiate the functions of these parts.

2. Calculate the total magnification produced by various combinations of eyepiece and objective lenses.
3. Demonstrate proper use of the microscope to observe and measure small objects.
4. Prepare a simple microscope slide and sketch the objects you observed.

EXPLORE

Procedure A—Microscope Basics

1. Familiarize yourself with the following list of rules for care of the microscope:
 - a. Keep the microscope under its *dustcover* and in a cabinet when it is not being used.
 - b. Handle the microscope with great care. It is an expensive and delicate instrument. To move it or carry it, hold it by its *arm* with one hand and support its *base* with the other hand (fig. 4.1).
 - c. Always store the microscope with the scanning or lowest-power objective in place. Always start with this objective when using the microscope.
 - d. To clean the lenses, rub them gently with *lens paper* or a high-quality cotton swab. If the lenses need additional cleaning, follow the directions in the "Lens-Cleaning Technique" section that follows.
 - e. If the microscope has a substage lamp, be sure the electric cord does not hang off the laboratory table where someone might trip over it. The bulb life can be extended if the lamp is cool before the microscope is moved.
 - f. Never drag the microscope across the laboratory table after you have placed it.
 - g. Never remove parts of the microscope or try to disassemble the eyepiece or objective lenses.
 - h. If your microscope is not functioning properly, report the problem to your laboratory instructor immediately.




Figure 4.1 Major parts of a compound light microscope with a monocular body and a mechanical stage. Some microscopes are equipped with a binocular body.

Lens-Cleaning Technique


1. Moisten one end of a high-quality cotton swab with one drop of Kodak lens cleaner. Keep the other end dry.
2. Clean the optical surface with the wet end. Dry it with the other end, using a circular motion.
3. Use a hand aspirator to remove lingering dust particles.
4. Start with the scanning objective and work upward in magnification, using a new cotton swab for each objective.
5. When cleaning the eyepiece, do not open the lens unless it is absolutely necessary.
6. Use alcohol for difficult cleaning and only as a last resort use xylene. Regular use of xylene will destroy lens coatings.

2. Observe a compound light microscope and study figure 4.1 to learn the names of its major parts. The lens system of a compound microscope includes three parts—the condenser, objective lens, and eyepiece or ocular.

Light enters this system from a *substage illuminator (lamp)* or *mirror* and is concentrated and focused by a *condenser* onto a microscope slide or specimen placed on the *stage*. The condenser, which contains a set of lenses, usually is kept in its highest position possible.

The *iris diaphragm*, located between the light source and the condenser, can be used to increase or decrease the intensity of the light entering the condenser. Locate the lever that operates the iris diaphragm beneath the stage and move it back and forth. Note how this movement causes the size of the opening in the diaphragm to change. (Some microscopes have a revolving plate called a disc diaphragm beneath the stage instead of an iris diaphragm. Disc diaphragms have different-sized holes to admit varying amounts of light.) Which way do you move the diaphragm to increase the light intensity? _____ Which way to decrease it? _____

After light passes through a specimen mounted on a microscope slide, it enters an *objective lens system*. This lens projects the light upward into the *body tube*, where it produces a magnified image of the object being viewed.



(a)



(b)

Figure 4.2 The powers of this 10× eyepiece (a) and this 40× objective (b) are marked in the metal. DIN is an international optical standard on quality optics. The 0.65 on the 40× objective is the numerical aperture, a measure of the light-gathering capabilities.

The *eyepiece (ocular) lens* system then magnifies this image to produce another image, seen by the eye. Typically, the eyepiece lens magnifies the image ten times (10×). Look for the number in the metal of the eyepiece that indicates its power (fig. 4.2). What is the eyepiece power of your microscope?

The objective lenses are mounted in a *rotating nosepiece* so that different magnifications can be achieved by rotating any one of several objective lenses into position above the specimen. Commonly, this set of lenses includes a scanning objective (4×), a low-power objective (10×), and a high-power objective, also called a high-dry-power objective (about 40×). Sometimes an oil immersion objective (about 100×) is present. Look for the number printed on each objective that indicates its power. What are the objective lens powers of your microscope?

To calculate the *total magnification* achieved when using a particular objective, multiply the power of the eyepiece by the power of the objective used. Thus, the 10× eyepiece and the 40× objective produce a total magnification of 10×40 , or 400×. See table 4.1.

Table 4.1 Microscope Lenses

Objective Lens Name	Common Objective Lens Magnification	Common Eyepiece Lens Magnification	Total Magnification
Scan	4×	10×	40×
Low power (LP)	10×	10×	100×
High power (HP)	40×	10×	400×
Oil immersion	100×	10×	1,000×

Note: If you wish to observe an object under LP, HP, or oil immersion, locate and then center and focus the object first under scan magnification.

- Complete Part A of Laboratory Report 4.
- Turn on the substage illuminator and look through the eyepiece. You will see a lighted circular area called the *field of view*.

You can measure the diameter of this field of view by focusing the lenses on the millimeter scale of a transparent plastic ruler. To do this, follow these steps:

- Place the ruler on the microscope stage in the spring clamp of a slide holder finger on a mechanical stage or under the stage (slide) clips. (*Note:* If

your microscope is equipped with a mechanical stage, it may be necessary to use a short section cut from a transparent plastic ruler. The section should be several millimeters long and can be mounted on a microscope slide for viewing.)

- b. Center the millimeter scale in the beam of light coming up through the condenser, and rotate the scanning objective into position.
- c. While you watch from the side to prevent the lens from touching anything, raise the stage until the objective is as close to the ruler as possible, using the *coarse adjustment knob* and then using the *fine adjustment knob* (fig. 4.3). (Note: The adjustment knobs on some microscopes move the body and objectives downward and upward for focusing.)
- d. Look into the eyepiece and use the coarse adjustment knob to raise the stage until the lines of the millimeter scale come into sharp focus.
- e. Adjust the light intensity by moving the *iris diaphragm lever* so that the field of view is brightly illuminated but comfortable to your eye. At the same time, take care not to overilluminate the field because transparent objects tend to disappear in bright light.
- f. Position the millimeter ruler so that its scale crosses the greatest diameter of the field of view.




Figure 4.3 When you focus using a particular objective, you can prevent it from touching the specimen by watching from the side.

Also, move the ruler so that one of the millimeter marks is against the edge of the field of view.

- g. In millimeters, measure the distance across the field of view.

5. Complete Part B of the laboratory report.
6. Most microscopes are designed to be *parfocal*. This means that when a specimen is in focus with a lower-power objective, it will be in focus (or nearly so) when a higher-power objective is rotated into position. Always center the specimen in the field of view before changing to higher objectives.

Rotate the low-power objective into position, and then look at the millimeter scale of the transparent plastic ruler. If you need to move the low-power objective to sharpen the focus, use the *fine adjustment knob*.

Adjust the iris diaphragm so that the field of view is properly illuminated. Once again, adjust the millimeter ruler so that the scale crosses the field of view through its greatest diameter, and position the ruler so that a millimeter mark is against one edge of the field.

Try to measure the distance across the field of view in millimeters.

7. Rotate the high-power objective into position while you watch from the side, and then observe the millimeter scale on the plastic ruler. All focusing using high-power magnification should be done only with the fine adjustment knob. If you use the coarse adjustment knob with the high-power objective, you can accidentally force the objective into the coverslip and break the slide. This is because the *working distance* (the distance from the objective lens to the slide on the stage) is much shorter when using higher magnifications.

Adjust the iris diaphragm for proper illumination. Usually more illumination when using higher magnifications will help you to view the objects more clearly. Try to measure the distance across the field of view in millimeters.

8. Locate the numeral 4 (or 9) on the plastic ruler and focus on it using the scanning objective. Note how the number appears in the field of view. Move the plastic ruler to the right and note which way the image moves. Slide the ruler away from you and again note how the image moves.
9. Observe a letter *e* slide in addition to the numerals on the plastic ruler. Note the orientation of the letter *e* using the scan, LP, and HP objectives. As you increase magnifications, note that the amount of the letter *e* shown is decreased. Move the slide to the left, and then away from you, and note the direction in which the observed image moves.
10. Examine the slide of the three colored threads using the low-power objective and then the high-power objective. Focus on the location where the three threads cross. By using the fine adjustment knob, determine the order from top to bottom by noting

which color is in focus at different depths. The other colored threads will still be visible, but they will be blurred. Be sure to notice whether the stage or the body tube moves up and down with the adjustment knobs of the microscope being used for this depth determination. The vertical depth of the specimen clearly in focus is called the *depth of field (focus)*. Whenever specimens are examined, continue to use the fine adjustment focusing knob to determine relative depths of structures clearly in focus within cells, giving a three-dimensional perspective. The depth of field is less at higher magnifications.



Critical Thinking Application

What was the sequence of the three colored threads from top to bottom? Explain how you came to that conclusion.

11. Complete Parts C and D of the laboratory report.

Demonstration

A compound light microscope is sometimes equipped with a micrometer scale mounted in the eyepiece. Such a scale is subdivided into fifty to one hundred equal divisions (fig. 4.4). These arbitrary divisions can be calibrated against the known divisions of a micrometer slide placed on the microscope stage. Once the values of the divisions are known, the length and width of a microscopic object can be measured by superimposing the scale over the magnified image of the object.

Observe the micrometer scale in the eyepiece of the demonstration microscope. Focus the low-power objective on the millimeter scale of a micrometer slide (or a plastic ruler) and measure the distance between the divisions on the micrometer scale in the eyepiece. What is the distance between the finest divisions of the scale in micrometers? _____




Figure 4.4 The divisions of a micrometer scale in an eyepiece can be calibrated against the known divisions of a micrometer slide. (Courtesy of Swift Instruments, Inc., San Jose, California)

Procedure B—Slide Preparation

1. Prepare several temporary *wet mounts* using any small, transparent objects of interest, and examine the specimens using the low-power objective and then a high-power objective to observe their details. To prepare a wet mount, follow these steps (fig. 4.5):
 - a. Obtain a precleaned microscope slide.
 - b. Place a tiny, thin piece of the specimen you want to observe in the center of the slide, and use a medicine dropper to put a drop of water over it. Consult with your instructor if a drop of stain might enhance the image of any cellular structures of your specimen. If the specimen is solid, you might want to tease some of it apart with dissecting needles. In any case, the specimen must be thin enough so that light can pass through it. Why is it necessary for the specimen to be so thin?
 - c. Cover the specimen with a coverslip. Try to avoid trapping bubbles of air beneath the coverslip by slowly lowering it at an angle into the drop of water.

- d. Remove any excess water from the edge of the coverslip with absorbent paper. If your microscope has an inclination joint, do not tilt the microscope while observing wet mounts because the fluid will flow.
 - e. Place the slide under the stage (slide) clips or in the slide holder on a mechanical stage, and position the slide so that the specimen is centered in the light beam passing up through the condenser.
 - f. Focus on the specimen using the scanning objective first. Next focus using the low-power objective, and then examine it with the high-power objective.
2. If an oil immersion objective is available, use it to examine the specimen. To use the oil immersion objective, follow these steps:
 - a. Center the object you want to study under the high-power field of view.
 - b. Rotate the high-power objective away from the microscope slide, place a small drop of immersion oil on the coverslip, and swing the oil immersion objective into position. To achieve sharp focus, use the fine adjustment knob only.
 - c. You will need to open the iris diaphragm more fully for proper illumination. More light is




Figure 4.5 Steps in the preparation of a wet mount.

needed because the oil immersion objective covers a very small lighted area of the microscope slide.

- d. The oil immersion objective must be very close to the coverslip to achieve sharp focus, so care must be taken to avoid breaking the coverslip or damaging the objective lens. For this reason, never lower the objective when you are looking into the eyepiece. Instead, always raise the objective to achieve focus, or prevent the objective from touching the coverslip by watching the microscope slide and coverslip from the side if the objective needs to be lowered. Usually when using the oil immersion objective, only the fine adjustment knob needs to be used for focusing.
3. When you have finished working with the microscope, remove the microscope slide from the stage and wipe any oil from the objective lens with lens paper or a high-quality cotton swab. Swing the scanning objective or the low-power objective into position. Wrap the electric cord around the base of the microscope and replace the dustcover.
4. Complete Part E of the laboratory report.

Demonstration

A stereomicroscope (dissecting microscope) (fig. 4.6) is useful for observing the details of relatively large, opaque specimens. Although this type of microscope achieves less magnification than a compound microscope, it has the advantage of producing a three-dimensional image rather than the flat, two-dimensional image of the compound light microscope. In addition, the image produced by the stereomicroscope is positioned in the same manner as the specimen, rather than being reversed and inverted as it is by the compound light microscope.

Observe the stereomicroscope. The eyepieces can be pushed apart or together to fit the distance between your eyes. Focus the microscope on the end of your finger. Which way does the image move when you move your finger to the right? _____

When you move it away? _____

If the instrument has more than one objective, change the magnification to a higher power. Use the instrument to examine various small, opaque objects available in the laboratory.




Figure 4.6 A stereomicroscope, also called a dissecting microscope.

11