





TRAINED STAFF MEMBERS

- 1. _____
- 2. _____
- 3. _____

- Room _____
- Room _____
- Room _____

<p>EpiPen® and EpiPen® Jr. Directions</p> <ul style="list-style-type: none">▪ Pull off gray activation cap.  <ul style="list-style-type: none">▪ Hold black tip near outer thigh (always apply to thigh).  <ul style="list-style-type: none">▪ Swing and jab firmly into outer thigh until Auto-Injector mechanism functions. Hold in place and count to 10. Remove the EpiPen® unit and massage the injection area for 10 seconds.	<p>Twinject™ 0.3 mg and Twinject™ 0.15 mg Directions</p>  <ul style="list-style-type: none">▪ Pull off green end cap, then red end cap.▪ Put gray cap against outer thigh, press down firmly until needle penetrates. Hold for 10 seconds, then remove. 
--	---

Once EpiPen or Twinject™ is used, call 911 or the Rescue Squad. Take the used unit with you to the Emergency Room. Plan to stay for observation at the Emergency Room for at least 4 hours.

For students with multiple food allergies, consider providing separate Action Plans for different foods.

**Medication checklist adapted from the Authorization of Emergency Treatment form developed by the Mount Sinai School of Medicine. Used with permission.